

Materials and Methods that Allow Fingerprint Analysis *and* DNA Profiling from the Same Latent Evidence


Step by Step Laboratory Protocol

Alexander Sinenikov (alex@ifi-test.com) and Karl A. Reich (karl@ifi-test.com)
Independent Forensics, 500 Waters Edge, Suite 210, Lombard, IL 60148, USA

Independent Forensics
DNA TESTING & TECHNOLOGIES

1. Collect *all* the Evidence

"Fingerprint & DNA" lift card




Lift card with hinged transparent liner. Easy opening tab on short side of card. Liner protects lifted print and allow easy access to biological material

Enhanced latent (starting material)



'Fingerprint' is lifted onto SIRCHIE/IFI lift card with transparent liner.




'Unlifted' biological material left behind on evidence.




'No DNA Left Behind'
Unlifted' biological material collected on swab. Retained for later extraction

2. Extract and Purify DNA from Adhesive Evidence (lift card, duct tape, sticky tape, etc.)

A) Immobilize lift card on bottom of lysis chamber



Remove liner from immobilizing circle



Place lift card on circle, adhesive layer facing up

B) Cover adhesive with Taming Material



Taming material being applied to adhesive



Adhesive layer is now hydrophilic


D) Combine all the biological fingerprint evidence



Use reserved swab from (1.), collect extract from Lysis Chamber

Assembled chamber is incubated at 56°C for 1 hour.

C) Extract DNA from Adhesive Layer



Add O-ring, cover and binder clips to form Lysis Chamber



Add lysis buffer to immobilized, covered fingerprint

E) Combined biological fingerprint evidence is extracted



Combined lysate and No DNA Left Behind material is recovered via spin-basket, and extracted at 56°C for 1 hour


F) DNA Purification on Xs (subtractive) Column



Note retained fingerprint powder on column, purified lysates

Can concentrate lysate (Centricon or SpeedVac), PCR amplify using kit of your choice.

GeneScan electropherogram of PowerPlex 16 PCR reaction after concentration of purified lysate.



3. Analyze ALL of your data; post PCR purification & concentration (AmpliconRx™)

PCR Reaction

+

Binding Buffer



load on column

Spin 3 min

[Discard run through]

Formamide + LIZ/ROX

[Incubate 5 min, RT °C]

Spin 2 min

Load Elution on CE

PCR samples that show partial profile are subjected to post-PCR purification and concentration using AmpliconRx™ kit

Up to 15X more CE signal

